

DoseASSURE™

OPTIMIZING INFLAMMATORY BOWEL DISEASE
TREATMENT THROUGH BIOLOGICS



DoseASSURE™, LabCorp's portfolio of biologics monitoring assays, may help physicians maximize treatment response using a personalized, patient-specific approach

- Help aid in titrating doses or adjusting frequency to optimize effectiveness¹⁻³
- May help avoid lack of response due to under-treatment¹
- Assist in preventing and managing loss of response due to immunogenicity⁴⁻⁵
- Minimize cost to patient by avoiding unhelpful dose escalation, especially in the setting of immunogenicity^{1,6}

Biologic Drug Name	Primary Target	*Clinical Indications	Test Name	Test No.
Infliximab Remicade®; Inflectra®, Renflexis®	TNF	CD, UC **	Infliximab and Anti-Infliximab Antibody (Serial Monitor), DoseASSURE™ IFX	503870
Adalimumab Humira®	TNF	CD, UC, RA	Adalimumab and Anti-Adalimumab Antibody (Serial Monitor), DoseASSURE™ ADL	503890
Vedolizumab Entyvio®	α4β7 integrin	CD, UC	Vedolizumab and Anti-Vedolizumab Antibody, DoseASSURE™ VDZ	504567
Golimumab Simponi®	TNF	UC, RA	Golimumab and Anti-Golimumab Antibody, DoseASSURE™ GOL	504563
Ustekinumab Stelara®	IL23, IL12	CD, PA, PP	Ustekinumab and Anti-Ustekinumab Antibody, DoseASSURE™ UST	504594
Certolizumab Cimzia®	TNF	CD, RA, PA, PP	Certolizumab and Anti-Certolizumab Antibody, DoseASSURE™ CTZ	504627

*Partial listing of FDA-approved indications. TNF: tumor necrosis factor, IL: interleukin, CD: Crohn's Disease Ulcerative Colitis, RA: Rheumatoid Arthritis, PA: Psoriatic Arthritis, PP: Plaque Psoriasis

**Also approved for pediatric forms of CD & UC

DoseASSURE test portfolio provides tests for both drug concentration (TDM) & anti-drug antibody (immunogenicity)

Therapeutic Drug Monitoring (TDM)

- Biologics have variable pharmacokinetics.^{3,7}
- Dosing by weight and empiric dose adjustment are inefficient and suboptimal.^{3,7}
- TDM for Biologics is a valuable tool to evaluate doses and to tailor dose adjustments to your individual patient.^{3,7}
- TDM can help differentiate under-treatment from other causes of lack of response.
- Proactive dose optimization using TDM may improve clinical scores and prolong duration of anti-TNF therapy.¹
- TDM has been shown to be cost-effective and may direct more appropriate care.¹⁶

Immunogenicity Testing (Anti-drug Antibody level)

- All biologics have the potential to induce an antibody-mediated immune response.
- Close to half of IBD patients on biologic therapy may develop anti-drug antibodies.^{4,8,9}
- Anti-drug antibodies may appear as early as after the first infusion and persist for years.⁸
- Anti-drug antibodies can adversely affect the amount of drug in the body.⁸
- Sufficient drug levels (e.g. infliximab >3ug/mL), concomitant use of immunomodulating agents, and regular dosing may protect against the risk of developing anti-drug antibodies.¹⁷⁻¹⁹

Interpreting Drug Concentrations

- Detectable drug levels are associated with better clinical outcome as measured by mucosal healing, lower C-reactive protein, higher remission rate, and less relapse.^{1,2,10,11}
- Target ranges and maximally effective concentrations have not been established.³
- Optimal drug concentration depends on the desired therapeutic endpoint and may differ case by case.¹²

Drug	Normal half-life	Proposed Target Trough Concentrations§
Infliximab	7.7 to 9.5 days	3 – 7 µg/mL ¹ ; 5 -10 µg/mL ² ; >4.0 µg/mL for mucosal healing ¹² ; ≥10.0 µg/mL may be required for fistula healing ²⁰
Adalimumab	Approx 2 weeks	≥ 7.5 µg/mL ¹³ > 5.85 µg/mL ¹⁴
Vedolizumab	Approx 25 days	>30 µg/mL at week 6 ¹³ >14 µg/mL during maintenance ²⁴
Golimumab	Approx 2 weeks	≥ 4.27 µg/mL correlated with greater response and remission ²²
Ustekinumab	Approx 3 weeks	>4.5 µg/mL has been associated with greater rate of endoscopic response ²³
Certolizumab	Approx 2 weeks	≥20 µg/mL correlated to higher remission rate ¹³

§Note: These target ranges were those used in landmark studies and do not necessarily translate into general recommendations for individual patients.

Interpreting Anti-Drug Antibody Levels

- Anti-drug antibodies can impact pharmacokinetics, efficacy, and the cost effectiveness of biologics.
- Low titer antibodies may have little to no effect on drug levels or clinical outcome but evidence suggests they may lead to later development of higher titers.⁹
- In contrast, high titers of antibodies are likely to be more consequential, leading to loss of drug efficacy by preventing drug binding to TNF and/or increasing drug clearance.^{9,15}
- Anti-drug antibody positivity should be interpreted in the context of the concomitant free drug level.

Anti-Drug Antibodies	Quantitative Range	Result Interpretation
Anti-Infliximab Abs	22- 10,000+ng/mL	Antibodies are reported as Low, Intermediate or High Titer
Anti-Adalimumab Abs	25-10,000+ ng/mL	Antibodies are reported as Low, Intermediate or High Titer
Anti-Vedolizumab Abs	25-10,000+ ng/mL	Stratification into low to high titer has yet to be determined.
Anti-Golimumab Abs	20-10,000+ ng/mL	Stratification into low to high titer has yet to be determined.
Anti-Ustekinumab Abs	40-10,000+ ng/mL	Stratification into low to high titer has yet to be determined.
Anti-Certolizumab Abs	40-10,000+ ng/mL	Stratification into low to high titer has yet to be determined.

- Patient-specific clinical context must be taken into account when evaluating drug and anti-drug antibody
- Serial measurements over time may be helpful

When & where to collect blood on my patients?

- The timing of sample collection is important because the drug concentration will change during the dosing interval.
- The Trough Concentration (TC) is measured at the least variable time in the dosing interval, just before the next dose (same day to within <7 days depending on the drug's normal half-life).
- During induction and maintenance phases, trough collections are usually recommended because target ranges are defined using TC.
- Blood can be drawn at any of LabCorp's nearly 2000 patient service centers located nationwide.

References

1. Vande Castele N, et al. Trough Concentrations of Infliximab Guide Dosing for Patients With Inflammatory Bowel Disease. *Gastroenter* 2015;148:1320-1329.
2. Vaughn BP, et al. Proactive Therapeutic Concentration Monitoring of Infliximab May Improve Outcomes for Patients with Inflammatory Bowel Disease: Results from a Pilot Observational Study. *Inflamm Bowel Dis* 2014;20:1996-2003.
3. Vaughn BP, et al. Biologic Concentration Testing in Inflammatory Bowel Disease. *Inflamm Bowel Dis* 2015;21:1435-4142.
4. Ungar B, et al. The temporal evolution of antidrug antibodies in patients with inflammatory bowel disease treated with infliximab. *Gut* 2014;63:1258-1264.
5. American Gastroenterological Association. Guidelines for the Identification, Assessment and Initial Medical Treatment in Crohn's Disease. <https://www.gastro.org/IBDcarepathway>
6. Steenholdt C, et al. Individualised therapy is more cost-effective than dose intensification in patients with Crohn's disease who lose response to anti-TNF treatment: a randomised, controlled trial. *Gut* 2014;63:919-927.
7. Ordas, et al. Therapeutic Drug Monitoring of Tumor Necrosis Factor Antagonists in Inflammatory Bowel Disease. *Clin Gastroenterol Hepatol* 2012;10:1079-1087.
8. Steenholdt C, et al. Clinical Implications of Variations in Anti-infliximab Antibody Levels in Patients with Inflammatory Bowel Disease. *Inflamm Bowel Dis* 2012 Vol 18(12):2209-2217
9. Vande Castele N, et al. Antibody Response to Infliximab and its Impact on Pharmacokinetics can be Transient. *Am J Gastroenterol* 2013;108:962-971.
10. Maser EA, et al. Association of Trough Serum Infliximab to Clinical Outcome After Scheduled Maintenance Treatment for Crohn's Disease. *Clin Gastroenterol Hepatol* 2006;4(10):1248-1254.
11. Seow CH, et al. Trough serum infliximab: a predictive factor of clinical outcome for infliximab treatment in acute ulcerative colitis. *Gut* 2010;59:49-54.
12. Imaeda H, et al. Relationship between serum infliximab trough levels and endoscopic activities in patients with Crohn's disease under scheduled maintenance treatment. *J Gastroenterol* 2014; 49:674-682.

13. Vande Castele N, et al. American Gastroenterological Association Institute Technical Review on the Role of Therapeutic Drug Monitoring in the Management of Inflammatory Bowel Diseases. *Gastroenterol*. 2017;153:835-857.
14. Mazar Y, et al. Adalimumab drug and antibody levels as predictors of clinical and laboratory response in patients with Crohn's disease. *Aliment Pharmacol Ther* 2014;40:620-628.
15. Yanai H, et al. Levels of Drug and Antidrug Antibodies Are Associated With Outcome of Interventions After Loss of Response to Infliximab or Adalimumab. *Clin Gastroenterol Hepatol* 2015;13(3):522-530.
16. Velayos F, et al. A Test-based Strategy Is More Cost-Effective than Empiric Dose Escalation for Patients with Crohn's Disease Who Lose Responsiveness to Infliximab 2013;11:654-666.
17. Colombel JR et al. Infliximab, Azathioprine, or Combination Therapy for Crohn's Disease. *N Engl J Med* 2010 362(15):1383-1395.
18. Brandse JR, et al. Insufficient Infliximab Exposure Predisposes to Immunogenicity and Enhanced Clearance of Infliximab in IBD. *J Crohns Colitis* 2016; 150(4 suppl 1): S144
19. Hanauer SB, et al. Incidence and Importance of Antibody Responses to Infliximab After Maintenance or Episodic Treatment in Crohn's Disease. *Clin Gastroenterol Hepatol* 2004; 2:542-553
20. Yarus A, et al. Higher Infliximab Trough Levels Are Associated With a Higher Rate of Perianal Fistula Healing in Patients With Crohn's Disease. *Gastroenterol* 2016;150(4):S105-S106.
21. Rosario M, et al. Relationship between vedolizumab pharmacokinetics and endoscopic outcomes in patients with UC. *Inflamm Bowel Dis* 2014;20:S1-53.
22. Sandborn WJ, et al. Subcutaneous Golimumab Induces Clinical Response and Remission in Patients With Moderate-to-Severe Ulcerative Colitis. *Gastroenterol* 2014;146:85-95.
23. Battat R, et al. Association Between Ustekinumab Trough Concentrations and Clinical, Biomarker, and Endoscopic Outcomes in Patients With Crohn's Disease. *Clin Gastroenterol Hepatol*. 2017;15:1427-1434.
24. Dreesen E, et al. Evidence to Support Monitoring of Vedolizumab Trough Concentrations in Patients With Inflammatory Bowel Diseases. *Clin Gastroenterol Hepatol*. 16(12):1937-46.



www.LabCorp.com

Visit www.LabCorp.com or call 800-444-9111 for more information.